Point of Care Certification Record Accu-Chek Inform II Glucose Meter

Name (First/Last)				Date		
Employee or Student #		Position (RN/ LPN/ Student)		Hospital / School		
COMPETENCY					ainer tials	Trainee Initials
Identifies base unit components						
Understands the color status						
Identifies meter components						
Checks battery status, communication status and date & time						
Understands annual recertification requirements:						
 4 QC (2 HI , 2 LO) and 5 Patient tests 						
Operator is locked out if certification expires						
Demonstrates QC testing						
Scans test strips & QC lot # correctly						
Aware to write opening date on QC vials; Discord 2 months often enoning						
 Discard 3 months after opening Inverts vial to mix, wipes tip & discards first drop 						
 Waits for flashing "drop" & applies sample drop to strip correctly 						
 Performs correct quality control steps with QC solution 						
Explains corrective action when QC FAILS & selects pre-defined comment						
Demonstrates patient testing						
Understands skin puncture process and supplies						
knows acceptable sample types Sinulates national (trins in matter hefere explains black)						
 Simulates patient testing/strip in meter before applying blood Understands to discard first drop of blood and use second drop for testing 						
 Explains critical values /actions /collection of lab samples for confirmation 						
 Knows to document results immediately on patient chart to prevent errors 						
Understands cleaning frequency & drying technique;						
Documents cleaning in meter						
Passed Accu-Chek Inform II quiz						
Knows location of on-line CBGM DSTs and training material (employees only)						
I understand the education presented today and will follow procedures and policies as outlined in the Point of Care Policy.						
Exam mark: (if < 85%, fill in action plan/ corrective action)						
Trainee (Signature): Date:						
Print as double-sided document Version: 1.1 Page 1 of 2						
Revision Date: 2014-0		version		Unique Identifi	er: I CBG	Page 1 of 2 0003 F1
Filename and Path: X:\QSE DOCUMENTS\LQSM-Finals\Point of Care\Accu-Chek Glucose Meter\I CBG 0003 F1 Point of Care Certification Record - Glucose Meter.docx						

MOST CURRENT VERSION ONLINE. UNCONTROLLED PRINTED COPY EXPIRES IN 24 HOURS

Trainer to complete this side of document

Certification Requirement:

- Trainee must be certified by a qualified TRAINER or EDUCATOR
- Trainee must complete written or on-line exam with a passing grade of >85% (18/20)
- This document must be signed by the trainee and a qualified educator

Qualified Educator:	Name (print)	
	Signature: Clinical Instructor/ Trainer/CPL	Date:

ACTION PLAN/ CORRECTIVE ACTION:

(Trainer to complete this section if 3 or more incorrect exam questions or any other competencies not met)

CERTIFICATION of TRAINERS / EDUCATORS

Extra training has been provided to this TRAINEE who is now qualified as a TRAINER or EDUCATOR to train others.

Job Title of new Educator: _____ DEPT: ____

Email a completed copy of this Certification Record and a copy of the exam to the laboratory.

Retained by POC Coordinator or Chief/Charge Technologist for documentation of certification.

 Print as double-sided document
 Version: 1.1
 Page 2 of 2

 Revision Date: 2014-08-25
 Unique Identifier:
 I CBG 0003 F1

 Filename and Path: X:\QSE DOCUMENTS\LQSM-Finals\Point of Care\Accu-Chek Glucose Meter\I CBG 0003 F1 Point of Care Certification Record - Glucose Meter.docx
 MOST CURRENT VERSION ONLINE. UNCONTROLLED PRINTED COPY EXPIRES IN 24 HOURS

YES